



Accord Medical Products Pvt. Ltd.

Manufacture, Exporter & Importers Of:
Physiotherapy, Occupational Therapy, Panchakarma Products,
Rehabilitation Aids, Hospital Furniture & Holloware, Orthopedic
Implants, Suction Machine & Slimming Equipments.

315, RG Complex-1 , Parshant Vihar, Rohini, Delhi-110085

Ph.No. - 09999303118, 09953888183

Off: 011-47543164

Email ID – info@accordmedicalproducts.com

Website – www.accordmedicalproducts.com


EXERCISE UNITS




AMP-03Ex76
Quadricep Table(withback rest)



AMP-03Ex77
Quadricep Table(without back rest)




AMP-03Ex62
Shoulder Wheel (Wall Mounting)



AMP-03Ex63
Axial Shoulder Wheel




AMP-03Ex71 Heel Exerciser



AMP-03Ex72 Ankle Exerciser

EXERCISE UNITS



AMP-03Ex65 T Pully




AMP-03Ex60 Shoulder Ladder



AMP-03Ex73 Rowing Machine



AMP-038082 Power Web(Size: 14")
AMP-038082 Power Web(Size: 7")



AMP-0326 Grip Exerciser (With Base)



AMP-03EX70 Leg Exerciser (Paddler)

EXERCISE UNITS



AMP:0325 Grip Exerciser
(Plastic with 4 Springs)



AMP:0320 Grip Exerciser
(V Shape)

AMP:0321 Grip Exerciser
(V Shape with Counter)



AMP:0330 Physio Ball (55cm, 76cm)

AMP:0334 Physio Ball (COSCO)
(55cm, 65cm, 75cm, 85cm,
95cm)

AMP-038094 Physio Ball (Gymnic)
(12",18",22",26",30",34",38",48")



AMP:038081 Gym Band (Imported)
(Red,Green,Blue,Silver,Golden)



AMP: 03EX20 Dumbbells
(1Kg to 5Kg)



AMP:0347 Weight Cuff
(0.5kg, 1kg, 1.5kg, 2kg, 2.5kg, 3kg)




AMP:0311 Foam Grip Ball

AMP:0314 Gel Grip Ball
(Small, Medium, Large)

AMP-0315 PVC Grip Ball
(small, medium, large)



AMP:032809 Shoulder Pulley Kit
(Over Door Model)



AMP -03EX20 Exercise Tubing
with handles (18",36",48")
light, medium, heavy, ultra heavy



AMP - 03EX61 Ball medicine
(1kg to 5kg)

AMP - 03EX62 Ball medicine
(Imported) 1/2kg to 5kg

Hot and Cold Therapy

AMP:03HC16 AIRCAST CRYO/CUFF COMPRESSION DRESSING SYSTEM

Consist of three parts:

- 1)A cuff that covers the specific body part with pressurized ice water.
- 2)A cooler that holds enough water and ice for 6 to 8 hours of cryotherapy.
- 3)A tube that exchanges the water between the cooler and cuff.

AutoChill system(cold separately,shown) eliminates manual water recycling by automatically exchanging water between cuff and cooler using an electric pump. The new IC Cooler features an integrated pneumatic pump with in the cooler lid that provides automated compression and cold therapy.simplicity of design and ease of operation makes it ideal for post-operative recovery,trauma,athletic training rooms and home use.



AMP:030816 Magic Warmer
(Imported, Size: 8"x16")



AMP:03HC04 Wax Bath
(Size: 18"x9"x7")



AMP:030507 Hot Cold Gel
Pack (Imported, Size: 5"x7")

AMP:030710 Hot Cold Gel
Pack (Imported, Size: 7"x10")



AMP:03HC09 Moist Heat
Therapy Unit (8Pack)



AMP:0345 Physio Wax (0.5 Kg)

AMP:0346 Physio Wax (1kg)



AMP:0346 Steam Pack

(Size: 5 x12, 10 x 12,
Contour, 10 x 22)



AMP:0351 Ultra Sound Gel (5kg)

AMP:0352 Ultra Sound Gel (250g)

Traction Aids



AMP-03TR01 Traction Unit



AMP-03TR02 Traction Bed(3 Fold)



AMP-03TR02 Traction Unit



AMP-03TR03 Traction Bed(4 Fold)

Treatment Furniture



AMP-03EX74 Tilt Table (Mannual)



AMP-03EX75 Tilt Table (Electrical)

Static cycle



AMP:03SL26: Static cycle(Economy)



AMP:03SL26A: Static cycle(Air Bike Commercial)



AMP:03SL26B Static cycle(Air Bike Semi.Commercial)



AMP:03SL26D: Static cycle(Single Shock)



AMP:03SL26C: Static cycle(Dual Shock)

AMP:03SL29 -Motorised Treadmill (with 6year Gaurantee on motor)



Motor Power : 3 HP DC Peak (Industrial Grade Motor)
Speed : 1 to 16 Km/hr
Incline : 0 to 15% (16 levels)
Running Belts : 20" * 58"
Programs : Preset, 10 levels, 1 Manual, 1HRC, 1
Backlit LCD with Speed/ Incline, Instant
Keys, Speed, Incline, Distance & Calories.

AMP:03T6300 -Motorised Treadmill (with 10year Gaurantee on motor)



Motor	:	3.0 HP AC
Drive System	:	Inverter Drive
Speed	:	1.0KM/HR-20KM/HR
Incline	:	0 – 20%
Running Area	:	520mm X 1530mm
Display	:	LED 6 windows
Pulse Sensor	:	Touch sensor
Max User Wt.	:	180 Kgs
Suspension	:	Smart suspended cushion
Deck	:	Elastic type
Belt	:	Nylon

AMP:03SL29A -Motorised Treadmill (Nexus) (with 6year Gaurantee on motor)



Motor Power : 3.0 HP AC (Continuous)
6.0 HP AC (Peak)
Speed : 0.8 to 22 Km/h
Incline : 0 to 15 % (16 Levels)
Running Belt : 22" x 60"
Programs : 10 Preset exercising programs with
5 Intensity Levels, 1 HRC, 2 Custom
& 4 User Programs. 3 LED &
1matrix, Time, Distance, Pulse,
Calories, Speed, Incline & Body Fat
Control System.

Active Passive Trainer

AMP:031267 Lower Extermity Active Passive Trainer



Sepcifications	
Set-up Size	130x60x105 cm
Motor	220V/DC
Frame	Steel
Shroud	ABS
Passive Mode	
Minimum	10 Rpm
Maximum	50 Rpm
Active Mode	
Resistance	6 levels

AMP:031266 Lower & Upper Extermity Motion Exercise Bike (Activ Passive Trainer)



DimensionsW	
Height	42
Width	60
Length	40-60
Package	56x25x105 cm
Passive Mode	
Minimum	10 Rpm
Maxmum	50 Rpm
Active Mode	
Resistance	6 levels
Certification	CE 93/42/EEC



AMP:031265 E-CYCLE ACTIVE/PASSIVE EXERCISE QUIPMENT

Description:

- ◌ ADVANTAGE: Body-Trainer to increase arm strength and maintain shoulder and elbow joint flexibility
- ◌ Foot-Trainer to improve endurance, strength, co-ordination and circulation
- ◌ SIZE:W51xD60xH40cm

AMP:03APT01 ACTIVE/PASSIVE Trainer (TZORA, ISRAEL)

Highly portable, weighing only 22 lbs / 10kg the APT-1 can be positioned anywhere. Place it on the floor for lower body exercises or on a table top for upper body workout.

Other features include exclusive, easily changeable accessories that enable complete utilization of the APT-1 by users with different and varied physical difficulties.

Passive Mode

5 speeds for passivE exercising. Auto reverse feature ensures release of seized muscles. Upper and lower limb training for passive, disabled users. Improve circulation, muscle tone and range of motion. Decrease atrophy, spasms and swelling.Reduces total limb resistance to passive cycling.Normalizes lower extremities muscle tone.

Active Mode

5 levels of resistance for Active exercising.Active training for upper and lower limbs.Adjustable crank-arm eight and length.Improve blood circulation, strength endurance and muscle-tone. Learn the motorics of cycling. Activate rhythmical locomotion patterns.





AMP-03SL24
Recumbent Bike



**AMP:03TE05: Compression
Therapy Unit**



AMP-03SL22
Crazy Disc
